

Retail Solutions

Volt CPAP Battery Pack



The portable Volt CPAP Battery Pack provides an alternative power supply and works with most major manufacturer's CPAP devices. It comes in a convenient compact carry case that makes it a portable and easy-to-carry power source — great for travel, camping and even power outages!

The Viverity Volt Battery Pack comes complete with battery, AC charger and handy carrying case. The carrying case has a handle, shoulder strap and four external zippered pouches, give it tremendous storage capability. The unit will provide power for up to 9 hours of operation at 15 cm H₂O of pressure and up to 14 hours at 10 cm H₂O. Designed with a universal 12-volt socket, it can accommodate standard cigarette lighter plugs from most manufacturer's power cords.

NOTE: The DC power cord for your CPAP device is not included with the Volt CPAP battery pack. DC power cords must be purchased specific to your brand. Please see the chart below or ask your representative for further information.




Product Specifications

- Comes complete with battery, AC charger and carrying case
- Provides power for up to 14 hours at 10 cm H₂O and 9 hours at 15 cm H₂O
- Product Weight: 2.5 lbs
- Charger: 15.6 VDC-3A output, 100-240 VAC input
- Battery Pack: 14.8 VAC
- Battery Type: Li-ion
- Rated Capacity: 7,800 mAh
- Minimum Capacity: 7,500 mAh
- Rated Energy: 115.44 volts
- Loaded Voltage: 14.4 volts
- Charge Voltage: 15.6 volts
- 90 day warranty



Recharging Power Cord



Carry Bag

Item #	Description	UPC Code	Qty
C100-4R	Volt CPAP Battery Pack	815067013203	each
VOLT-BAG	Replacement Carry Bag for the Volt CPAP Battery	092237604915	each



To Order:

(800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com