

## CPAP Accessories

# RemZzzs CPAP Mask Liners

The RemZzzs CPAP mask liner provides a soft layer between the patient's face and the silicone found on most CPAP mask cushions. The liner can help alleviate several common CPAP therapy side effects. RemZzzs CPAP liners are made from a special blend of natural fibers which help prevent leaks and irritation around the cushion seal while absorbing facial oils and moisture.



## Product Specifications

- Created to fit most popular CPAP masks
  - Made from a special blend of 100% naturally absorbent cotton fibers
  - Acts as a barrier between skin and the silicone of the mask, helping to eliminate most of the common problems associated with wearing a CPAP mask
  - Significantly reduces air leaks, allowing more consistent pressure
  - Prevents skin irritation and ugly pressure marks
  - Helps to prevent the cushion from slipping by absorbing facial oils which tend to loosen the mask seal
  - Assist with comfortably holding patient's mask in place
  - Promotes compliance
  - One fresh facial liner per night is recommended for maximum effectiveness and comfort
  - A forehead liner is provided, for use with masks with a forehead cushion
  - Beards and mustaches do not interfere with the effectiveness of RemZzzs
  - In order to offer the highest comfort, RemZzzs will extend about a 1/2 to 3/4 of an inch beyond mask cushion's edge
  - Each 30-day kit comes with 30 RemZzzs mask liners and 2 forehead liners
- NOTE: each forehead liner lasts for about 15 days

Step 1



Place the RemZzzs directly on top of the mask cushion.

Step 2



Lean downward towards mask, placing nose and mouth into oval shaped hole.

Step 3



Press face against RemZzzs while returning to upright position.

Step 4



Snap mask fasteners into place and adjust tightness.

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:

(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com

## Mask Specifications

### Full Face Masks

Item #	Size	Mask	Qty
K1-FS	Small	Roscoe Full Face Mask, small Roscoe DreamEasy™ Full Face Mask, small Roscoe Sapphire Full Face Mask, small Resmed Quattro/FX, small Resmed Ultra Mirage, small Philips Respironics ComfortGel, small Philips Respironics ComfortFull2, small Philips Respironics FullLife, small EVO Comfort-Fit, APEX	30 liners per box, 2 forehead pads
K2-FM	Medium	Roscoe Full Face Mask, medium Roscoe DreamEasy™ Full Face Mask, medium Roscoe Sapphire Full Face Mask, medium Resmed Quattro/FX, medium and large Resmed Ultra Mirage, medium Philips Respironics ComfortGel, medium Philips Respironics ComfortFull2, medium Philips Respironics FullLife, medium EVO Comfort-Fit, APEX	30 liners per box, 2 forehead pads
K3-FL	Large	Roscoe Full Face Mask, large Roscoe DreamEasy™ Full Face Mask, large Roscoe Sapphire Full Face Mask, large Resmed Quattro/FX, large Resmed Ultra Mirage, large Philips Respironics ComfortGel, large Philips Respironics ComfortFull2, large Philips Respironics FullLife, large EVO Comfort-Fit, APEX	30 liners per box, 2 forehead pads
K4-FS	Small	Fisher & Paykel FlexiFit 431, small Fisher & Paykel FlexiFit 432, small Fisher & Paykel Forma, small Hans Rudolph 7600 V2, small	30 liners per box, 2 forehead pads
K5-FM	Medium	Fisher & Paykel FlexiFit 431, medium Fisher & Paykel FlexiFit 432, medium Fisher & Paykel Forma, medium Hans Rudolph 7600 V2, medium	30 liners per box, 2 forehead pads
K6-FL	Large	Fisher & Paykel FlexiFit 431, large Fisher & Paykel FlexiFit 432, large Fisher & Paykel Forma, large Hans Rudolph 7600 V2, large	30 liners per box, 2 forehead pads
6C-AVK	One Size	Amara View	30 liners per box

### Nasal Masks

Item #	Size	Mask	Qty
K7-NS	Small	Philips Respironics ComfortGel, small Philips Respironics ComfortClassic, small Philips Respironics Profile Lite, small Philips Respironics Comfort Select, small Philips Respironics Comfort Fusion, small	30 liners per box, 2 forehead pads
K8-NM	Medium	Philips Respironics ComfortGel, medium Philips Respironics ComfortClassic, medium Philips Respironics Profile Lite, medium Philips Respironics Comfort Select, medium Philips Respironics Comfort Fusion, medium	30 liners per box, 2 forehead pads
K9-NL	Large	Philips Respironics ComfortGel, large Philips Respironics ComfortClassic, large Philips Respironics Profile Lite, large Philips Respironics Comfort Select, large	30 liners per box, 2 forehead pads
K10-NS	Small	Roscoe DreamEasy™, small Roscoe Nasal Mask, small Roscoe Sapphire Nasal Mask, small Resmed Activa and Activa LT, small Resmed Mirage Micro, small Resmed Ultra Mirage II Shallow/Wide, small Fisher & Paykel FlexiFit 405, small Fisher & Paykel FlexiFit 406, small Fisher & Paykel Zest Petite, small	30 liners per box, 2 forehead pads
K11-NM	Medium	Roscoe DreamEasy™, medium Roscoe Nasal Mask, medium Roscoe Sapphire Nasal Mask, medium Resmed Activa and Activa LT, medium Resmed Mirage Micro, medium Resmed Ultra Mirage II Standard, medium Fisher & Paykel FlexiFit 405, medium Fisher & Paykel FlexiFit 406, medium Fisher & Paykel Zest One Size, medium	30 liners per box, 2 forehead pads
K12-NL	Large	Roscoe DreamEasy™, large Roscoe Nasal Mask, large Roscoe Sapphire Nasal Mask, large Resmed Activa and Activa LT, large Resmed Mirage Micro, large Resmed Ultra Mirage II, large Fisher & Paykel Zest Plus, large	30 liners per box, 2 forehead pads

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:

(800) 3-ROSCOE  
fax: (440) 572-4261  
roscoemedical.com