



Aluminum Crutches

ProBasics Aluminum Crutches are sturdy, lightweight anodized aluminum crutches with push-button height adjustment. The extra thick underarm pads and grips offer comfort for use. The double extruded center tube is a solid one-piece design providing additional strength to weight-bearing area. Non-skid rubber tips provide excellent traction and absorb shock. Easy wing nut adjustments for hand grips. Available in youth, adult or tall sizes to accommodate all individuals. Supports patient weights up to 350 pounds.



Extra thick underarm pads



Easy push-button height adjustment

FEATURES

- Made of lightweight, anodized aluminum
- Reinforced one-piece center tube provides additional strength to weight-bearing area
- Push-button adjustment makes changing height easy
- Height adjustment in 1-inch increments
- Hand grip position readily adjusts with wing nuts
- Non-skid rubber tips provide excellent traction on many surfaces
- Extra thick underarm pads and grips offer additional comfort
- HCPCS Code: E0114*

SPECIFICATIONS

Patient Height Range:	Junior: 4'6" – 5'2" Adult: 5'2" – 5'10" Tall: 5'10" – 6'6"
Arm Pad to Hand Grip:	Junior: 9" – 13" Adult: 15" – 19" Tall: 16" – 20"
Arm Pad to Floor:	Junior: 37" – 46" Adult: 45" – 53" Tall: 53" – 61.5"
Weight Capacity:	350 lb
Product Weight:	Junior: 3.3 lb Adult: 3.6 lb Tall: 3.8 lb
Carton Dimensions:	Junior: 39" x 13" x 11" Adult: 47" x 13" x 11" Tall: 54" x 13" x 11"
Shipping Weight:	Junior: 29 lb Adult: 32 lb Tall: 34 lb
Warranty:	Limited lifetime on frame 90-day on all parts

Item #	Description	UPC	UOM
CRJ	ProBasics Aluminum Crutches, Junior	815067071289	8/cs
CRA	ProBasics Aluminum Crutches, Adult	815067071272	8/cs
CRAT	ProBasics Aluminum Crutches, Tall	815067071296	8/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.