



Shower Chair

The **ProBasics Shower Chair** is designed for individuals who require seating support and stability when bathing or showering. The chair features a comfortable contoured textured seat and back that can support patient weights up to 250 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Blow-molded bench and back provide comfort and strength
- Quick tool-free assembly
- Wide-based crutch tips help increase stability
- Two-hole handles in seat provide security when sitting or standing
- Legs are height adjustable with one-inch snap button adjustments
- HCPCS Code: E0240*

SPECIFICATIONS

Overall Footprint:	16.5"W x 19"D
Seat Dimensions:	19.5"W x 11.5"D
Adjustable Seat Height:	14.25" – 20.5"
Weight Capacity:	250 lb
Carton Dimensions:	28" x 21" x 19"
Shipping Weight:	25 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSCWB	ProBasics Shower Chair with Back	815067071913	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.